

# Application Note A1

## IviumStat and TDI-load operation



### ***Impedance measurements on a 3-cell sealed lead acid battery using an IviumStat and an electronic load***

#### **1. Basic principles**

A potentiostat/galvanostat is used to control and measure electrochemical cells. The potentiostat/galvanostat is capable of "linear control" of the electrochemical cell from negative to positive potentials and currents (and vice versa), i.e. charging and discharging of a battery. However, most potentiostat/galvanostats are limited in the power/current output. In some of these cases it is possible to use a programmable electronic load to test high power electrochemical systems. An electronic load in general has much higher current capacity.

In principle a load is similar, but not as versatile, compared to a potentiostat/ galvanostat: loads can only "consume" power and not deliver it, and loads generally do not have a frequency response analyser available. Moreover, loads, by operating principle, cannot go "through zero". In the battery example: loads cannot charge a battery, nor can loads short circuit a battery down to 0V potential. But in various cases a load is useful for testing purposes when high currents are required, i.e. for battery and fuel cell testing.

To accommodate testing with an electronic load and an available potentiostat, Ivium potentiostats can be operated to control the load, and do frequency response measurements on an electrochemical cell with/via the load.

#### **2. Device description**

An Ivium potentiostat/galvanostat, like the IviumStat can control peripheral equipment in two ways, i.e. by using the peripheral outputs or using the cell output. Because all Ivium potentiostats/galvanostats have a built-in frequency response analyser, also impedance measurements can be done in combination with a load. Our standard IviumStat rates  $\pm 5A @ \pm 10V$ , with various potential and current ranges.

A load can vary in shapes, sizes and capacities. Depending on the make and model they come with different input/output interfaces. For this experiment a TDI-Dynaload was used, model RBL488-800, with a maximum voltage of 100 V (ranges 10-50-100V full scale) and a maximum current of 120 A (2-20-120A full scale). It was equipped with an IEEE-488 interface (RS-232

optional but not installed on this device) and an analog interface (terminal block).

In principle, higher capacity loads and different make&model loads are operated in a similar manner.

As test subject a sealed lead acid (LA) battery was used (HYS613). The LA battery consisted of 3 cells with a nominal voltage of 6 V and a nominal capacity of 1.3 Ah.

#### **3. Connection**

The TDI RBL488 is equipped with a 9-points terminal block:

##### Terminal Block on RBL488

CS	→	Current sample output 0-10V, 0-Fullscale (CS tov S-)
REM	→	Remote program input External prog. Source. 0-10V, 0-Fullscale (REM tov S-) <sup>1</sup>
SYN	→	Sync output
EN	→	Remote enable Connection between EN-S- will turn on the DC
S-	→	Negative Sense
S-	→	Negative Sense
E-	→	Negative Input
S+	→	Positive Input
E+	→	Positive Input

The connections between the IviumStat and the Dynaload were made as shown below:

<u>IviumStat cell cable</u>		<u>RBL488 terminal</u>
CE+RE <sup>2</sup>	→	REM
S+GND	→	S
WE+10kOhm <sup>3</sup>	→	CS

The positive and negative leads of the Dynaload were connected respectively to the positive and negative terminals of the sealed LA battery.

<sup>1</sup>The externally applied signal is deducted from 10V, the result is applied to the cell.

The optimal way to use the capability of the IviumStat and IviumSoft software is using the cell cable directly to control the load.

<sup>2</sup>Because the CE+RE is negative w.r.t. the set potential at the WE, a negative potential needs to be set in order to get a positive potential at the REM.

<sup>3</sup>A 10 kOhm resistor was connected between WE and CS to accommodate the maximum output current of the CS terminal. With a 10 V maximum output of CS and the 10 kOhm resistor in place the maximum current measured on WE is 1 mA, which is taken as the current range for the impedance measurement.




Figure 1: Schematic of the connections IviumStat-Dynaload

## 4. Experimental

With the equipment set-up as described measurements were carried out on the LA battery. Because a load can only consume power the battery needed to be operated in the discharge region, so a bias voltage was applied.

### 4.1 LSV

First an LSV was recorded with the IviumStat controlling the Dynaload, then the IviumStat was connected directly to the LA battery (via the standard cell cable). The experiments were carried out in CV mode. Because of the control of the Dynaload the values chosen in the IviumSoft need to be different;

for the Dynaload operation:

- E-Start = -3.95 V, E-end = -4.5 V, will result in a 6.05 V to 5.5 V scan
- E step = 5 mV, Scanrate = 0.025 V/s
- IviumStat current range CR=1mA
- 4-Electrode cell connection (default)

for the IviumStat operation:

- E-Start = 6.05 V, E-end = 5.5 V,
- E step = 5 mV, Scanrate = 0.025 V/s
- IviumStat current range CR=1A, will be measured directly
- 4-Electrode cell connection (default)

### 4.2 Impedance

An impedance experiment was carried out with the IviumStat controlling the Dynaload, then the IviumStat was connected directly to the LA battery and the experiment was repeated. The method parameters for each measurement are given in Annex 1. As can be seen from these parameters for the Dynaload operation the values need to be chosen differently;

for the Dynaload operation:

- DC-bias: -4 V, to get +4V at the CE+RE and 10-4=6V will be applied to the cell
- IviumStat current range CR=1mA
- 4-Electrode cell connection (default)

for the IviumStat operation:

- DC-bias: 6V, will be applied to the cell directly
- IviumStat current range CR=1A, will be measured directly

- 4-Electrode cell connection (default)

In both measurements the frequency was scanned in single sine mode from 5000Hz down to 3Hz with 10 frequencies per decade and an amplitude of 20mV (parameters: see Annex 1). The Dynaload was operated in CV mode.

## 5. Results

### 5.1 LSV

The results for the LSV scan are given in Figure 2 <sup>4</sup>.




Figure 2: LSV results

As can be seen in these results, both measurement directly with an IviumStat and measurement using an IviumStat to control the Dynaload yield stable results. Note that the difference in current between the experiments is caused by the difference in state of charge of the LA battery.

### 5.2 Impedance

The results for the impedance measurements are given in the graphs below <sup>5</sup>.

<sup>4</sup> When comparing the LSV results, a multiplying factor needs to be taken into account: given that the load was operated at 20A full scale, and full scale equals an output of 10V, which will result in a measurement of 1 mA at the WE at full scale, this factor equals 2E+4.

<sup>5</sup> When comparing the impedance results, a multiplying factor needs to be taken into account: given that the load was operated at 120A full scale, and full scale equals an output of 10V, which will result in a measurement of 1 mA at the WE at full scale, this factor equals 1.2E+5.




Figure 3: Impedance results

As shown in the graph above, both the measurement via the load and the measurement directly with the IviumStat show the same results up to ca. 1 kHz frequency. It should be noted that the battery was operated at a DC bias in the discharge region, leading to a bias discharge current. This explains the variance in the impedance as a function of frequency (and time).

At higher frequencies, the load shows odd behaviour, the IviumStat continues to perform well. This can be explained by looking at the phase angle during the measurement (see Figure 4).




Figure 4: Phase angle results

In figure 4 it shows that the phase angle in the frequency region 3 Hz to 1E+3 Hz is about equal for both devices. However, when using the load it shows that towards higher frequencies the absolute phase angle increases. This is expected, given that the bandwidth of the load is lower (ca. 1.6 kHz; determined from the slew rate setting of 100µs) than that of the IviumStat (1 MHz).

These results show that in this particular case, the load should not be operated at frequencies over ~1 kHz <sup>6</sup>.

### Annex 1: Method parameters impedance measurement

	Value	Unit
Title	Scan 1	
E start	-4.0000	V
Equilibration time	0	s
Frequencies	33	
Current Range	1mA	
+Noise Reduction	<input type="checkbox"/> Off	
Filter	automatic	
Stability	high speed	
Connect to	Cell-4EL	
Analog inputs	none	
+AutoCR	<input type="checkbox"/> Off	
+DualCR	<input type="checkbox"/> Off	
+Apply wrt DCP	<input type="checkbox"/> Off	
+Cell after meas	<input type="checkbox"/> Off	
Pretreatment	0	levels
Data Options	<input type="checkbox"/>	
+Automatic save	<input type="checkbox"/> Off	
+AUX	<input type="checkbox"/> Off	
+Anout2	<input type="checkbox"/> Off	
MeasConfig	standard	
Report	<input type="checkbox"/>	

Method parameters for measurement via load

	Value	Unit
Title	Scan 1	
E start	6.0000	V
Equilibration time	0	s
Frequencies	33	
Current Range	1A	
+Noise Reduction	<input type="checkbox"/> Off	
Filter	10kHz	
Stability	high speed	
Connect to	Cell-4EL	
Analog inputs	none	
+AutoCR	<input type="checkbox"/> Off	
+DualCR	<input type="checkbox"/> Off	
+Apply wrt DCP	<input type="checkbox"/> Off	
+Cell after meas	<input type="checkbox"/> Off	
Pretreatment	0	levels
Data Options	<input type="checkbox"/>	
+Automatic save	<input type="checkbox"/> Off	
+AUX	<input type="checkbox"/> Off	
+Anout2	<input type="checkbox"/> Off	
MeasConfig	standard	
Report	<input type="checkbox"/>	

Method parameters for measurement via IviumStat

<sup>6</sup> Note that different test objects may allow the use of different method parameters. In this particular case the rather small LA battery limited the use of higher currents, which otherwise might have yielded more stable test results.

Method frequencies:

The 'Edit frequencies' dialog box is shown with the 'Multi sine' tab selected. The 'Frequencies' section is set to '10' with an 'Apply' button. The 'Frequency scan' section has 'Start' at 5000 Hz, 'End' at 3 Hz, and 'Amplitude' at 0.02 V. A 'Manual override' checkbox is unchecked. Below these settings is a table with 14 rows of frequency data.

	freq / Hz	amp / mV
freq[1]	5000.00	20.00
freq[2]	3965.40	20.00
freq[3]	3144.88	20.00
freq[4]	2494.14	20.00
freq[5]	1978.06	20.00
freq[6]	1568.76	20.00
freq[7]	1244.15	20.00
freq[8]	986.71	20.00
freq[9]	782.54	20.00
freq[10]	620.62	20.00
freq[11]	492.20	20.00
freq[12]	390.35	20.00
freq[13]	309.58	20.00
freq[14]	245.52	20.00