# **Application Note A5**

## **RRDE** rotator



### Introduction

Ivium is currently developing a rotator for a rotating ring disk (RRDE) electrode, which can, for example, be used to study reaction mechanisms and ionic diffusion coefficients related to redox chemistry. The rotator has dual brush contacts, the rotation rates are variable between 100 and 10,000 RPM in steps of 1 RPM and it is compatible with Ivium tips as well as 3<sup>rd</sup> party tips. The rotation speed can be controlled manually from the manual control unit, or via analog output of Ivium potentiostats.

Recently, Ivium performed measurements with the RRDE rotator. In this report the results are presented.

#### Experimental

A CompactStat.h standard was used in combination with the RRDE rotator. The rotation speed was set manually with the standalone control unit. The measurements were carried out in a threeelectrode configuration. The working electrode was an RRDE platinum ring/glassy carbon disk electrode (4 mm) from ALS, Japan. The counter electrode (CE) and reference electrode (RE) were a platinum sheet and an Ag/AgCl electrode (3M NaCl), respectively. The rotator, counter electrode and reference electrode were positioned in a 100 ml glass beaker with a standard laboratory stand. The measurements were carried out in a 0.1 M KNO3 aqueous solution containing 10 mM K<sub>3</sub>Fe(CN)<sub>6</sub> at room temperature. A picture of the experimental setup is shown in Figure 1.



Figure 1: Photograph of the RRDE setup.

Cyclic voltammetry measurements with a scan rate of 100 mV/s and a step size of 10 mV were performed between -1.0 and 1.0 vs Ag / AgCl at various rotation rates.

#### **Results and discussion**

Figure 2 shows the results of the current response as a function of electrode rotation rate. As expected, the reduction current increases by increasing the rotation rate.



Figure 2: Cyclic voltammograms at various rotation rates at a scan rate of 100 mV/s and a step size of 10 mV.

The results obtained in the cyclic voltammetry measurements can be used to determine the diffusion coefficient of the  $Fe(CN)_6^{3-}$  ion by using the Levich equation, which is given by

$$i_l = 0.62 n FAD^{2/3} \omega^{1/2} v^{-1/6} C_0$$

where  $i_l$  is the limiting current under mass transfer controlled conditions, n is the number of electrons in the redox reaction, F is the Faraday constant, Dis the diffusion coefficient,  $\omega$  is the rotation rate, vis the kinematic viscosity of the electrolyte and  $C_0$  is the concentration of the ion in the bulk of the solution.<sup>1</sup>

Figure 3 shows that the current density is indeed proportional to square root of the rotation rate as expected based on the Levich equation. From the slope of the linear fit the diffusion coefficient of the  $Fe(CN)_6^{3-}$  ion is calculated to be  $5.92 \times 10^{-6}$  cm<sup>2</sup>/s. This value is close to known values such as  $7.26 \times 10^{-6}$  cm<sup>2</sup>/s in 1 M KCl at  $25^{\circ}$ C.<sup>2</sup> Therefore, it can be concluded that the rotator works well.



Figure 3: Levich graph of the limiting current as a function of rotation rate.

#### Conclusions

Preliminary results of the RRDE rotator show that the limiting current is proportional to the rotation rate satisfies the Levich equation. The calculated diffusion coefficient of the  $Fe(CN)_6^{3-}$  ion is close to reported values.

#### References

- 1. A.J. Bard and L.R. Faulkner, *Electrochemical methods, fundamentals and applications*, 2nd Ed., John Wiley & Sons, Inc., 2001.
- 2. S.J. Konopka and B. McDuffie, *Diffusion* coefficients of ferri- and ferrocyanide ions in aqueous media, using twin-electrode thin-layer electrochemistry, Anal. Chem., 42 (1970) 1741-1746.